
Ion Exchange Chromatography





Sepax Technologies, Inc.

Sepax Technologies, Inc. develops and manufactures products in the area of chemical and biological separations, bio-surfaces and proteomics. Sepax product portfolio includes 1) liquid chromatography columns and media, 2) SPE and Flash chromatography columns and tubes, 3) bulk resin for preparative separation and process chromatography, and 4) natural product and Chinese traditional medicine separation and purification.



A leader in Biological Separations

Sepax develops and manufactures wide range of biological separation products using both silica and polymeric resins as the support. The selection of particle size is from 1 µm to 100 µm and pore size from non-porous to 2000 Å. Unique and proprietary resin synthesis and surface technologies have been developed for solving the separation challenges in biological area.



Size Exclusion

 $SRT^{\mathbb{B}}$

SRT®-C

Nanofilm[®]

ZenixTM

ZenixTM-C

Ion-exchange

Proteomix[®]

Glycomix TM

Antibody Separation

AntibodixTM

Carbohydrate Separation

Carbomix[®]

Analytical, Semi-prep and Preparative



Proteomix[®] Ion-exchange Phases

General Description

Proteomix ion-exchange columns are specially designed for high resolution, high efficiency and high recovery separations of proteins, oligonucleotides, carbohydrates, and peptides. The packing support is composed of nonporous, rigid, spherical, highly cross-linked poly (styrene divinylbenzene) (PS/DVB) beads. The non-porous resins have particle sizes of 1.7, 3, 5 and 10 μm . The PS/DVB resin surface is grafted with a highly hydrophilic, neutral polymer layer with the thickness in the range of nanometer. The hydrophobic PS/DVB resin surface is totally covered by such a hydrophilic coating that eliminates non-specific bindings with biological analytes, leading to high efficiency and high recovery separations for bio-molecules. On the top of the hydrophilic layer, a layer of ion-exchange functional groups is attached via chemical bonding. A proprietary chemistry was developed to synthesize a densely packed, uniform ion-exchange layer.

Chemical Structure of Proteomix Resins

The chemical structure of Proteomix ion-exchange phases is composed of a rigid PS/DVB core, a densely packed, nanometer thick, hydrophilic coating, and a uniform ion-exchange layer, as shown in Figure 1.




Figure 1. Schematic illustration of the chemical structure of Proteomix ion-exchange phases.

Proteomix Ion-Exchange Phases

As shown in Figure 2, Proteomix ion-exchange phases include SCX, WCX, SAX, and WAX, which are strong cation exchanger with sulfonate functional groups, weak cation exchanger with carboxylate functional groups, strong anion exchanger with quaternary ammonium functional groups, weak anion exchanger with tertiary amine functional, respectively. All ion-exchange functional groups are chemically bonded to the top of the hydrophilic coating.




Figure 2. Chemical compositions of Proteomix SCX, WCX, SAX, and WAX phases.

Highlights of Proteomix Non-Porous Ionexchange Resins

- Combined high capacity and high separation efficiency in one material
- Particle size selection of 1.7, 3, 5 and 10 μm
- Mono-dispersed particles
- \bullet Tolerate high pressure operation: 4,000, 6,000, 8,000 and 12,000 psi for 10, 5, 3 and 1.7 μm resins, respectively
- Wide pH range: 2-12
- High resolving power for slightly differed structures of biological species
- \bullet 1.7 and 3 μm particles are best suitable for whole protein separation in proteome studies
- Suitable for separations of peptides, carbohydrates and glycans, proteins, polynucleotides, cell lysates and multidimensional separations


Non-Porous Proteomix Ion-exchangers – Technology Advancements

High Capacity of Proteomix Non-Porous Resins

It is well known that the non-porous resin maximizes the mass transfer and minimizes the lateral diffusion, resulting in high speed and high efficiency separations. However, due to the low surface area, the separation capacity is low. The low capacity issue creates a great challenge to the non-porous resins for many applications. The Proteomix non-porous ion-exchange resins developed by Sepax proprietary technologies have a breakthrough on the improvement of their capacities. This new technology enables the capacity of the non-porous resins increased to the level comparable to that of the porous resins. The dynamic binding capacity data are listed in the technical specifications.

High binding capacity leads to high sample loading without deterioration of the separation efficiency and the resolution. Figure 3 shows the elution of cytochrome C with the loading up to 200 μ g for a 1.7 μ m Proteomix SCX non-porous column without any deterioration of the separation.

Figure 3. Elution profiles of cytochrome C at various loadings.



Column: Proteomix SCX-NP1.7 (1.7 µm, 4.6x50 mm)

Mobile Phase: A, 20 mM phosphate buffer, pH 6.0;

B, A+1.0 M NaCl

Gradient: 0-70%B (21 min)
Flow Rate: 0.35 mL/min
Detection: 280 nm
Temperature: 25 °C


Sample: Cytochrome C (20 mg/mL)

High Separation Efficiency, Resolution and Selectivity

Proteomix SCX-NP, WCX-NP, SAX-NP, and WAX-NP resins have three unique features. First, their nanometer thick hydrophilic layer completely eliminates the non-specific interactions with the biological analytes. Second, non-porous beads minimize the biological analytes' lateral

diffusion and suppress their diffusion into the chromatographic bed. Third, a uniform layer of ion-exchange functional groups synthesized by Sepax's proprietary technology greatly improves their ion-exchange capacities. Such innovatively designed Proteomix SCX-NP, WCX-NP, SAX-NP, and WAX-NP phases result in the highest efficiency separations for proteins, oligonucleotides, carbohydrates, and peptides. Figure 4 is a typical test chromatogram for separation of three proteins: ribonuclease A, cytochrome C, and lysozyme by a 4.6×50 mm, Proteomix SCX-NP column (3 μ m). The efficiency of lysozyme reaches 100,000 of plates with 5 cm long column. Such a high efficiency separation is unprecedented.

Figure 4. Separation of a protein mixture by a Proteomix SCX-NP column.



Column: Proteomix SCX-NP3 (3 µm, 4.6x50 mm)

Mobile phase: A, 10 mM phosphate, pH 6.0

B, A + 1.0 M NaCl

Gradient: 0-70%B in 15 min

Flow rate: 0.5 mL/min

Sample: 1) Cytochrome C, 2) Ribonuclese A,

3) Lysozyme

Injection: $5 \mu L (1 \text{ mg/mL for each protein})$

Temperature: 25 °C Detection: UV 280 nm

The uniqueness of non-porous Proteomix SCX, WCX, SAX, and WAX phases offers highest resolution and selectivity for protein separations. Figure 5 shows an elution profile of lysozyme by a Proteomix SCX-NP3 phase (3 μm). Six peaks are well resolved with a short column (4.6 \times 50 mm). Such a high selectivity and resolving power is unmatched with any other SCX phases.

Figure 5. The elution profile of lysozyme from its multiple impurities by a Proteomix SCX-NP column (3 μ m, 4.6 \times 50 mm). The separation conditions are the same as those in Figure 6.





Figure 6 is a typical test chromatogram for a 5 μ m, Proteomix SAX-NP column for separation of a mixture of ovalbumin and BSA. High resolution and high selectivity of Proteomix SAX-NP phase can well separate the impurities contained in ovalbumin sample, as well as the BSA dimer from BSA.

Figure 6. Elution profile of a mixture of ovalbumin and BSA by Proteomix SAX-NP5 column.



Column: Proteomix SAX-NP5 (5 µm, 4.6x50 mm)

Mobile phase: A, 20 mM Tris, pH 8.0

B, A + 0.5 M NaCl

Gradient: 0-50%B in 15 min

Flow rate: 0.5 mL/min

Sample: 1) Ovalbumin (10 mg/mL), 2) BSA (5


mg/mL)

Injection: $10 \mu L$ Temperature: $25 \,^{\circ}C$ Detection: UV 214 nm

Charge Dependent Separation

Figure 7 is a typical test chromatogram for a Proteomix SCX-NP column (3 μ m). Extremely high efficiency separation was achieved for 5 proteins. Those proteins were eluted based on their charge status (pI) from low to high: ovalbumin, ribonuclease A, cytochrome C, aprotinin, and lysozyme. At the separation conditions (pH 6.0), ovalbumin is negatively charged and eluted at the void volume, indicating no non-specific bindings.

Figure 7. Separation of a protein mixture with various pI by a Proteomix SCX column.



Column: Proteomix SCX-NP3 (3 µm, 4.6x50 mm)

Mobile phase: A, 10 mM phosphate, pH 6.0

B, A + 1.0 M NaCl

Gradient: 0-100%B in 15 min

Flow rate: 0.5 mL/min


Injection: 5 µL (1 mg/mL for each protein)

Temperature: 25 °C Detection: UV 280 nm

High Lot-to-Lot Reproducibility

With well-controlled resin production and the surface chemistry, the manufacturing of the Proteomix ion-exchangers is highly reproducible. The typical variation of the retention time is less than 6% from batch to batch. One example is shown in Figure 8 for the production of three lots of 3 μ m Proteomix SAX-NP resins.

Figure 8. The reproducibility of three lots of Proteomix SAX resins.



Column: Proteomix SAX-NP3 (3 µm, 4.6x50 mm)

Mobile phase: A, 20 mM Tris, pH 8.0

B, A + 0.3 M NaCl


Gradient: 0-100%B in 15 min

Flow rate: 0.5 mL/min

Sample: 1) Ovalbumin (10 mg/mL), 2) BSA (5

mg/mL)

Injection: $10 \mu L$ Temperature: $25 \,^{\circ}C$ Detection: UV 214 nm



Mono-dispersed Proteomix ionexchange particles enable uniform column packing bed, which improves separation efficiency.

Column Dimension Availability

The column dimensions of the Proteomix SCX, WCX, SAX, and WAX products are 0.75, 2.1, 3.0, 4.6, 7.8, 10, and 21.2 mm I.D., and 2, 3, 5, 10, 15, 25, and 30 cm length. The columns are available both in stainless steel and PEEK tubing. Sepax also offers custom-made columns.



Features

- ♦ Uniform polymer beads as the packing support
- Proprietary surface chemistry specially designed for elimination of non-specific bindings
- Unprecedented separation efficiency, selectivity and resolving power
- Complete selection for analytical, semi-preparative and preparative separations
- ♦ Wide particle size selection
- ♦ High stability
- ♦ High recovery
- Well suited for UPLC system (1.7 μm)
- Compatible with FPLC system (5 μm and 10 μm)

Technical Specifications

Products	Pore Size	Particle Size (µm)	Dynamic binding capacity*	pH range
Proteomix SCX-NP1.7, NP3, NP5 & NP10	Non-porous	1.7, 3, 5, 10	~30, 22, 13, 10 mg/mL	2-12
Proteomix WCX-NP1.7, NP3, NP5 & NP10	Non-porous	1.7, 3, 5, 10	~21, 18, 16, 7.4 mg/mL	2-12
Proteomix SAX-NP1.7, NP3, NP5 & NP10	Non-porous	1.7, 3, 5, 10	~34, 31, 17, 13 mg/mL	2-12
Proteomix WAX-NP1.7, NP3, NP5 & NP10	Non-porous	1.7, 3, 5, 10	~25, 23, 16, 14 mg/mL	2-12

^{*}Dynamic binding capacity test conditions: 1) For Proteomix SCX and Proteomix WCX, Sample: 3.0 mg/mL lysozyme in 20 mM sodium phosphate buffer, pH 6.0; Flow rate: 0.5 mL/min (0.25 mL/min for 1.7 μm); Detection: UV 254 nm. 2) For Proteomix SAX and Proteomix WAX, Sample: 3.0 mg/mL BSA in 20 mM Tris HCl, pH 8.0; Flow rate: 0.5 mL/min (0.25 mL/min for 1.7 μm); Detection: UV 280 nm .

Applications


Separation and Analysis					
Proteins					
Cell lysates					
Nucleic acids					
Peptides					
Carbohydrates and glycans					
Nucleotides					

Proteomix ion-exchange phases have wide applications for separation, identification and purification of proteins, protein variants, peptide fragments, phosphorylated, sialylated, deaminated, PEGylated and other derivatized proteins. They are excellent for monitoring enzymatic reactions and protein-protein interactions. With extremely high separation efficiency, Proteomix phases are well designed for proteomics applications and 2D or multi-dimensional separations. One example is the high resolution of cell lysates by Proteomix SAX-NP phases. Proteomix SAX and WAX phases are suitable for separation of oligonucleotides.

Separation of Horse Serum, a Protein Mixture

Figure 9 shows the separation of a horse serum by using Proteomix SAX-NP5 and a commercial porous SAX column. Non-porous Proteomix SAX-NP5 has much high resolution than the porous SAX column.

Figure 9. The separation of a horse serum (20 μ L, 2x diluted) from BioWhittaker, a Cambrex company (Walkersville, MD). (Courtesy of Miyako Kawakatsu, M&S Instruments, Inc.)



Mobile phase: (A) 25 mM Tris acetate, pH 7.3

(B) A + 0.8 M sodium acetate

Gradient: 0-50%-70% B (0-25-28 min)
Flow rate: 0.5 mL/min for SAX-NP5


1.0 mL/min for Porous SAX

Detection: UV 280 nm

Separation of Carbohydrates and Glycans

Figure 10 shows the separation of glucose oligomers and corn starch hydrolysate. ELSD is utilized to detect carbohydrate molecules. For ELSD detection, bleeding is a big problem for most columns, especially silica based amino columns. However, Proteomix SAX columns not only solved the bleeding problem, but also achieved excellent separation efficiency and resolution. Figure 11 is another example of high resolution separation of glycans and their isomers by non-porous Proteomix SAX column.

Figure 10. Separation of glucose oligomers and corn starch hydrolysate by Proteomix SAX-NP5 phase. (Courtesy of Miyako Kawakatsu, M&S Instruments, Inc.)



Column: Proteomix SAX-NP5 (5 μm, 4.6x150 mm)

Mobile phase: A, 0.05% (25% NH₄OH) in CH₃CN

B, 0.05% (25% NH₄OH) in H₂O


Gradient: 15-80% B (20 min)

Flow rate: 0.7 mL/min

Sample: Corn starch hydrolysate

Temperature: 35 °C Detection: ELSD

Figure 11. Separation of 2-AA (anthranilic acid) labeled N-linked oligosaccharides profiling of an IgG1 sample.



Column: Proteomix SAX-NP5 (5 μ m, 4.6x150 mm) Mobile phase: A, 2.5% (v/v) acetic acid, 0.5% TEA in H₂O

B, 0.5% acetic acid in ACN

Gradient: 0-100% B (60 min)
Flow rate: 0.3 mL/min


Detection: Fluorescence Ex/Em=360/425nm

Sample: G0F: asialo, agalacto, core-fucosylated biantennary glycan; G1F: asialo, mono-galacto, core-fucosylated biantennary glycan; G2F: asialo, di-galacto, core-fucosylated biantennary glycan; SA1F: mono-sialylated, galactosylated, core-fucosylated biantennary glycan; SA2F: di-sialylated, galactosylated, core-fucolylated biantennary glycan; NGNA: N-glycolylneuraminic acid.

Separation of Peptides in Volatile Buffer

The peptides with different hydrophobicity can be resolved by reversed phases, such as C18. However, peptides with different charges may not be well separated on reversed phase mode if their hydrophobic properties are close. As an alternative method, ion exchange chromatography is recommended. With the net charge from +1 to +4, four peptides (C1, C2, C3 and C4 with their sequences listed in Figure 12) were well separated by Proteomix SCX-NP phase in a volatile buffer of ammonium acetate and acetonitrile, which makes it compatible for LC/MS application.

Figure 12. Cation exchange separation of 4 peptides C1, C2, C3 and C4 (0.1 mg/mL). (Courtesy of Miyako Kawakatsu, M&S Instruments, Inc.)



Column: Proteomix SCX-NP3 (3 μm, 4.6x50 mm)

Mobile phase: A, 5 mM CH₃COONH₄:CH₃CN=4:1 (v/v)

B, 500 mM CH₃COONH₄:CH₃CN=4:1(v/v)

Gradient: 0-50%B (20 min)
Flow rate: 0.6 mL/min
Temperature: 25 °C
Detection: ELSD

Injection: $5 \mu L (0.1 \text{ mg/mL})$


Sample: Peptides C1, C2, C3 and C4

PeptideSequenceNet ChargeC1Ac-Gly-Gly-Gly-Leu-Gly-Gly-Ala-Gly-Gly-Leu-Lys-amide+1C2Ac-Lys-Tyr-Gly-Leu-Gly-Gly-Ala-Gly-Gly-Leu-Lys-amide+2C3Ac-Gly-Gly-Ala-Leu-Lys-Ala-Leu-Lys-Gly-Leu-Lys-amide+3C4Ac-Lys-Tyr-Ala-Leu-Lys-Ala-Leu-Lys-Gly-Leu-Lys-amide+4

High Protein Loading Separation

Ribonuclease A was separated by a non-porous SCX column. When the injection amount of Ribonuclease A was increased from 50 µg to 5 mg, the separation resolution had negligible change and each of the impurities was well separated, as shown in Figure 13. The combination of high loading capacity and high resolution separation enables non-porous Proteomix ion-exchange columns very well suited for both analytical and preparative separations of biological molecules.

Figure 13. Elution profiles of RNaseA by Proteomix SCX-NP3 column with various loadings. (Courtesy of Miyako Kawakatsu, M&S Instruments, Inc.)



Column: Proteomix SCX-NP3 (3 µm, 4.6x50 mm)

Eluent: 10 mM PBS, pH 6.0 Gradient: 0-1.0 M NaCl (20 min)

Flow rate: 1.0 mL/min
Detection: 280 nm;
Temperature: 25 °C
Injection: 100 µL


Sample: Ribonuclease A (from Bovine Pancreas for

Genomics, Wako Pure Chemical industries, Ltd.)

Column Screening - Particle Size Selection

Proteomix ion-exchange resins have a wide range of the particle size selection from 1.0 to 10 μm for analytical and semi-preparative separation applications. Figure 14 is the separation of three proteins by four weak cation exchange resins with the particle size of 1.7, 3, 5, and 10 μm , separately. It is clearly seen that the smaller particle generates higher efficiency and resolution. One of the impurities in lysozyme labeled by Peak N is barely separated with10 μm WCX resin. For 5 μm WCX resin, it is almost a baseline separation. When the particle size is decreased to 3 μm , it becomes a baseline separation. With the particle size further decreased to 1.7 μm , the impurity is well separated.

Figure 14. Separation of a protein mixture by Proteomix WCX-NP columns with various particle sizes.



Column: Proteomix WCX-NP (4.6x50 mm)
Mobile phase: A, 10 mM phosphate, pH 6.0

B, A + 1.0 M NaCl

Gradient: 0-100% B in 20 min

Flow rate: 1.0 mL/min, 0.75 mL/min for WCX-NP1.7

Sample: 1) Cytochrome C, 2) Ribonuclese A,

3) Lysozyme


Injection: $5 \mu L (1 \text{ mg/mL for each protein})$

Temperature: 25 °C Detection: UV 280 nm

Separation of Oligonucleotides

High resolution separation is critical for DNA analysis and purification. The non-porous Proteomix SAX resins offer high separation efficiency with high capacity. Figure 15 shows the separation of a 50-mer oligonucleotide synthesized by ABI synthesizer. The largest peak is the modified 50-mer oligo and the others are failures or the impurities. As a comparison, the same sample was separated with a commercial porous SAX column. The results showed similar retention time for the 50-mer oligonucleotide from both non-porous and porous SAX columns, which indicated that their capacities are close to each other. However, the non-porous SAX generated much higher efficiency and resolution.

Figure 15. Separation a 50-mer oligonucleotide from its impurities by non-porous Proteomix SAX and porous SAX columns.



Column: $4.6x150 \text{ mm } (5 \text{ } \mu\text{m})$


Mobile phase: A, 25 mM Tris, pH 9; B, A + 1.0 M LiCl

Gradient: 30-75% B in 45 min

Flow rate: 1.0 mL/min Sample: 50-mer Oligo Injection: 50 μ L Detection: UV 260 nm

The advantage of Proteomix SAX column in impurity analysis of oligonucleotides is again demonstrated in Figure 16. Compared with Vendor D anion-exchange column, Proteomix SAX column fully resolves impurities in the 21-mer oligonucleotide sample.

Figure 16. Separation of a 21-mer oligonucleotide: comparison of Proteomix SAX and Vendor D column.



Panel A conditions

Column: Proteomix SAX-NP5 (5 μ m, 4.6x250 mm) Mobile phase: A: 20 mM Tris (pH 8.0); B: A+0.5 M NaCl

Gradient: 0-70% B in 21 min Flow rate: 0.5 mL/min Wavelength: 260 nm Column temperature: Ambient

Panel B conditions

Column: Vendor D anion-exchange column (4.0x250 mm)


Mobile phase: A: 10mM NaClO₄ 1mM Tris pH 7.4

B: 300mM NaClO₄ 1mM Tris pH 7.4

Gradient: 0-75% B in 30 min Flow rate: 1 mL/min Column temperature: 50 °C

Figure 17 is another example of high resolution separation of oligonucleotides of Mw 12196 from its degraded fragments. At least 35 species of those degraded oligo fragments were well resolved within 22 minutes.

Figure 17. Separation of an oligonucleotide (MW 12196) and its degraded fragments.



Column: Proteomix SAX-NP5 (5 µm, 4.6x150 mm)
Mobile phase: A, 25 mM Tris, 1.0 mM EDTA, 10% ACN,

pH 8.0; B, A+1.0 M NaCl

Gradient: 0-75%B in 50 min Flow rate: 0.5 mL/min


Sample: Oligonucleotide and its degraded fragments


Injection: 50 µL

Separation and Analysis of Cell Lysates

The key issue for proteome studies is to separate and identify a large number of biological species in a cell, such as proteins, nucleotides, peptides and others. The demand for separation is unprecedented. With the uniqueness of high resolution and high capacity, the non-porous Proteomix ion-exchange resins are very much suitable for separating cell lysates. Figure 18a showed the separation profiles of *E. coli* lysate with 3, 5, and 10 μ m non-porous Proteomix SAX particles. The minimum number of resolved peaks increased from 40 to 60 to 75 when the particle size decreased from 10 to 5 to 3 μ m. To better view the separation performance, the elution profiles in the range of 10-22 minutes were shown in Figure 18b. At least 45, 38, 27 peaks were resolved at the retention time of 10-22 minutes for 3, 5, and 10 μ m SAX columns, respectively.

Figure 18. (a, Top) Separation of *E. coli* lysate by 3, 5 and 10 μ m Proteomix SAX particles. (b, Bottom) Expanded region of retention time from 10 to 22 min.





Column: Proteomix SAX-NP (4.6x50 mm)

Mobile phase: A, 20 mM Tris, pH 9.0; B, A + 0.5 M NaCl

Gradient: 0-100% B in 30 min

Flow rate: 0.5 mL/min

Sample: E. coli lysate (2.5 mg/mL)

Injection: $10 \mu L$ Temperature: $25 \,^{\circ}C$ Detection: UV 280 nm Figure 19 shows various sample loading for a 3 μ m, 4.6x50 mm non-porous SAX column. When the amount of *E. coli* lysate increased from 25 μ g to 50 μ g to 125 μ g, the separation efficiency and resolution remained consistent.

Figure 19. Elution profile of *E. coli* Lysate by non-porous Proteomix SAX column with various loadings. The separation conditions are the same as those in Fig. 18.





Figure 20. High peak capacity separation of yeast lysate by non-porous Proteomix SAX column. (Courtesy of Prof. Fred Regnier, Purdue University)



Column: Proteomix SAX-NP3 (3 µm, 4.6x50 mm)

Mobile phase: A, 20 mM Tris, pH 9.0; B, A + 0.5 M NaCl

Gradient: 0-100% B in 30 min

Flow rate: 0.5 mL/min

Sample: Yeast lysate (2.5 mg/mL)

Injection: $40 \,\mu L$ Temperature: $25 \,^{\circ}C$ Detection: UV 280 nm

Product Information of Proteomix Ion-Exchange Columns

Proteomix SCX-NP Columns

Phase	ID x Length (mm)	P/N	Phase	ID x Length (mm)	P/N
	7.8 x 50	401NP2-7805		7.8 x 50	401NP3-7805
	4.6 x 150	401NP2-4615		4.6 x 150	401NP3-4615
	4.6 x 100	401NP2-4610		4.6 x 100	401NP3-4610
	4.6 x 50	401NP2-4605		4.6 x 50	401NP3-4605
	4.6 x 30	401NP2-4603		4.6 x 30	401NP3-4603
	4.0 x 10 (Guard)	401NP2-4001	2	4.0 x 10 (Guard)	401NP3-4001
1.7 μm	3.0 x 100	401NP2-3010	3 μm Proteomix	3.0 x 100	401NP3-3010
Proteomix	3.0 x 50	401NP2-3005	SCX-NP3	3.0 x 50	401NP3-3005
SCX-NP1.7	3.0 x 30	401NP2-3003	30X-W 3	3.0 x 30	401NP3-3003
	2.1 x 100	401NP2-2110		2.1 x 100	401NP3-2110
	2.1 x 50	401NP2-2105		2.1 x 50	401NP3-2105
	2.1 x 30	401NP2-2103		2.1 x 30	401NP3-2103
	2.0 x 10 (Guard)	401NP2-2001		2.0 x 10 (Guard)	401NP3-2001
	Precolumn Filter (PEEK)*	102000-P356		Precolumn Filter (PEEK)	102000-P356
	10 x 250	401NP5-10025		10 x 250	401NP10-10025
	7.8 x 250	401NP5-7825		7.8 x 250	401NP10-7825
	7.8 x 100	401NP5-7810		7.8 x 150	401NP10-7815
	7.8 x 50	401NP5-7805		7.8 x 50	401NP10-7805
	4.6 x 250	401NP5-4625		4.6 x 250	401NP10-4625
	4.6 x 150	401NP5-4615		4.6 x 150	401NP10-4615
	4.6 x 100	401NP5-4610		4.6 x 100	401NP10-4610
	4.6 x 50	401NP5-4605		4.6 x 50	401NP10-4605
5 μm	4.6 x 30	401NP5-4603	10 μm	4.6 x 30	401NP10-4603
Proteomix	4.0 x 10 (Guard)	401NP5-4001	Proteomix SCX-NP10	4.0 x 10 (Guard)	401NP10-4001
SCX-NP5	2.1 x 150	401NP5-2115		2.1 x 150	401NP10-2115
	2.1 x 100	401NP5-2110		2.1 x 100	401NP10-2110
	2.1 x 50	401NP5-2105		2.1 x 50	401NP10-2105
	2.1 x 30	401NP5-2103		2.1 x 30	401NP10-2103
	2.0 x 10 (Guard)	401NP5-2001		2.0 x 10 (Guard)	401NP10-2001
	Precolumn Filter (PEEK)	102000-P355		Precolumn Filter (PEEK)	102000-P355
	Preparative C	Columns		Preparative Columns	
	21.2 x 250	401NP5-21225		21.2 x 250	401NP10-21225
	21.2 x 150	401NP5-21215		21.2 x 150	401NP10-21215

Proteomix WCX-NP Columns

Phase	ID x Length (mm)	P/N	Phase	ID x Length (mm)	P/N
	7.8 x 50	402NP2-7805	3 μm Proteomix WCX-NP3	7.8 x 50	402NP3-7805
	4.6 x 150	402NP2-4615		4.6 x 150	402NP3-4615
	4.6 x 100	402NP2-4610		4.6 x 100	402NP3-4610
	4.6 x 50	402NP2-4605		4.6 x 50	402NP3-4605
	4.6 x 30	402NP2-4603		4.6 x 30	402NP3-4603
4.7	4.0 x 10 (Guard)	402NP2-4001		4.0 x 10 (Guard)	402NP3-4001
1.7 μm Proteomix	3.0 x 100	402NP2-3010		3.0 x 100	402NP3-3010
WCX-NP1.7	3.0 x 50	402NP2-3005		3.0 x 50	402NP3-3005
WCX-NP1.7	3.0 x 30	402NP2-3003		3.0 x 30	402NP3-3003
	2.1 x 100	402NP2-2110		2.1 x 100	402NP3-2110
	2.1 x 50	402NP2-2105		2.1 x 50	402NP3-2105
	2.1 x 30	402NP2-2103		2.1 x 30	402NP3-2103
	2.0 x 10 (Guard)	402NP2-2001		2.0 x 10 (Guard)	402NP3-2001
	Precolumn Filter (PEEK)	102000-P356		Precolumn Filter (PEEK)	102000-P356

Phase	ID x Length ((mm)	P/N	Phase	ID x Length (mm)	P/N
	10 x 250	402NP5-10025		10 x 250	402NP10-10025
	7.8 x 250	402NP5-7825		7.8 x 250	402NP10-7825
	7.8 x 100	402NP5-7810		7.8 x 150	402NP10-7815
	7.8 x 50	402NP5-7805		7.8 x 50	402NP10-7805
	4.6 x 250	402NP5-4625		4.6 x 250	402NP10-4625
	4.6 x 150	402NP5-4615		4.6 x 150	402NP10-4615
	4.6 x 100	402NP5-4610		4.6 x 100	402NP10-4610
	4.6 x 50	402NP5-4605		4.6 x 50	402NP10-4605
5 μm	4.6 x 30	402NP5-4603	10 μm	4.6 x 30	402NP10-4603
Proteomix	4.0 x 10 (Guard)	402NP5-4001	Proteomix	4.0 x 10 (Guard)	402NP10-4001
WCX-NP5	2.1 x 150	402NP5-2115	WCX-NP10	2.1 x 150	402NP10-2115
	2.1 x 100	402NP5-2110		2.1 x 100	402NP10-2110
	2.1 x 50	402NP5-2105		2.1 x 50	402NP10-2105
	2.1 x 30	402NP5-2103		2.1 x 30	402NP10-2103
	2.0 x 10 (Guard)	402NP5-2001		2.0 x 10 (Guard)	402NP10-2001
	Precolumn Filter (PEEK)	102000-P355		Precolumn Filter (PEEK)	102000-P355
	Preparative Columns			Preparative	Columns
	250 × 21.2	402NP5-21225		21.2 x 250	402NP10-21225
	150 × 21.2	402NP5-21215		21.2 x 150	402NP10-21215

Proteomix SAX-NP Columns

Phase	ID x Length (mm)	P/N	Phase	ID x Length (mm)	P/N
	7.8 x 50	403NP2-7805		7.8 x 50	403NP3-7805
	4.6 x 150	403NP2-4615		4.6 x 150	403NP3-4615
	4.6 x 100	403NP2-4610		4.6 x 100	403NP3-4610
	4.6 x 50	403NP2-4605		4.6 x 50	403NP3-4605
	4.6 x 30	403NP2-4603		4.6 x 30	403NP3-4603
	4.0 x 10 (Guard)	403NP2-4001		4.0 x 10 (Guard)	403NP3-4001
1.7 μm	3.0 x 100	403NP2-3010	3 μm	3.0 x 100	403NP3-3010
Proteomix SAX-NP1.7	3.0 x 50	403NP2-3005	Proteomix SAX-NP3	3.0 x 50	403NP3-3005
070X141 1.7	3.0 x 30	403NP2-3003	0/0/1410	3.0 x 30	403NP3-3003
	2.1 x 100	403NP2-2110		2.1 x 100	403NP3-2110
	2.1 x 50	403NP2-2105		2.1 x 50	403NP3-2105
	2.1 x 30	403NP2-2103		2.1 x 30	403NP3-2103
	2.0 x 10 (Guard)	403NP2-2001		2.0 x 10 (Guard)	403NP3-2001
	Precolumn Filter (PEEK)	102000-P356		Precolumn Filter (PEEK)	102000-P356
	10 x 250	403NP5-10025		10 x 250	403NP10-10025
	7.8 x 250	403NP5-7825		7.8 x 250	403NP10-7825
	7.8 x 100	403NP5-7810		7.8 x 150	403NP10-7815
	7.8 x 50	403NP5-7805		7.8 x 50	403NP10-7805
	4.6 x 250	403NP5-4625		4.6 x 250	403NP10-4625
	4.6 x 150	403NP5-4615		4.6 x 150	403NP10-4615
	4.6 x 100	403NP5-4610		4.6 x 100	403NP10-4610
	4.6 x 50	403NP5-4605		4.6 x 50	403NP10-4605
5 μm	4.6 x 30	403NP5-4603	10 μm	4.6 x 30	403NP10-4603
Proteomix	4.0 x 10 (Guard)	403NP5-4001	Proteomix SAX-NP10	4.0 x 10 (Guard)	403NP10-4001
SAX-NP5	2.1 x 150	403NP5-2115		2.1 x 150	403NP10-2115
	2.1 x 100	403NP5-2110		2.1 x 100	403NP10-2110
	2.1 x 50	403NP5-2105		2.1 x 50	403NP10-2105
	2.1 x 30	403NP5-2103		2.1 x 30	403NP10-2103
	2.0 x 10 (Guard)	403NP5-2001		2.0 x 10 (Guard)	403NP10-2001
	Precolumn Filter (PEEK)	102000-P355		Precolumn Filter (PEEK)	102000-P355
	Preparative C			Preparative	
	21.2 x 250	403NP5-21225		21.2 x 250	403NP10-21225
	21.2 x 150	403NP5-21215		21.2 x 150	403NP10-21215

Proteomix WAX-NP Columns

Phase	ID x Length (mm)	P/N	Phase	ID x Length (mm)	P/N
	7.8 x 50	404NP2-7805		7.8 x 50	404NP3-7805
	4.6 x 150	404NP2-4615		4.6 x 150	404NP3-4615
	4.6 x 100	404NP2-4610		4.6 x 100	404NP3-4610
	4.6 x 50	404NP2-4605		4.6 x 50	404NP3-4605
	4.6 x 30	404NP2-4603		4.6 x 30	404NP3-4603
	4.0 x 10 (Guard)	404NP2-4001		4.0 x 10 (Guard)	404NP3-4001
1.7 μm	3.0 x 100	404NP2-3010	3 μm	3.0 x 100	404NP3-3010
Proteomix	3.0 x 50	404NP2-3005	Proteomix	3.0 x 50	404NP3-3005
WAX-NP1.7	3.0 x 30	40NP2-3003	WAX-NP3	3.0 x 30	404NP3-3003
	2.1 x 100	404NP2-2110		2.1 x 100	404NP3-2110
	2.1 x 50	404NP2-2105		2.1 x 50	404NP3-2105
	2.1 x 30	404NP2-2103		2.1 x 30	404NP3-2103
	2.0 x 10 (Guard)	404NP2-2001		2.0 x 10 (Guard)	404NP3-2001
	Precolumn Filter (PEEK)	102000-P356		Precolumn Filter (PEEK)	102000-P356
	10 x 250	404NP5-10025		10 x 250	404NP10-10025
	7.8 x 250	404NP5-7825		7.8 x 250	404NP10-7825
	7.8 x 100	404NP5-7810		7.8 x 150	404NP10-7815
	7.8 x 50	404NP5-7805		7.8 x 50	404NP10-7805
	4.6 x 250	404NP5-4625		4.6 x 250	404NP10-4625
	4.6 x 150	404NP5-4615		4.6 x 150	404NP10-4615
	4.6 x 100	404NP5-4610		4.6 x 100	404NP10-4610
	4.6 x 50	404NP5-4605		4.6 x 50	404NP10-4605
5 μm	4.6 x 30	404NP5-4603	10 μm	4.6 x 30	404NP10-4603
Proteomix	4.0 x 10 (Guard)	404NP5-4001	Proteomix WAX-NP10	4.0 x 10 (Guard)	404NP10-4001
WAX-NP5	2.1 x 150	404NP5-2115		2.1 x 150	404NP10-2115
	2.1 x 100	404NP5-2110		2.1 x 100	404NP10-2110
	2.1 x 50	404NP5-2105		2.1 x 50	404NP10-2105
	2.1 x 30	404NP5-2103		2.1 x 30	404NP10-2103
	2.0 x 10 (Guard)	404NP5-2001		2.0 x 10 (Guard)	404NP10-2001
	Precolumn Filter (PEEK)	102000-P355		Precolumn Filter (PEEK)	102000-P355
	Preparative C			Preparative	
	21.2 x 250	404NP5-21225		21.2 x 250	404NP10-21225
	21.2 x 150	404NP5-21215		21.2 x 150	404NP10-21215

^{*} Precolumn Filters comes with 0.5 μ m PEEK frit for 102000-P356 and 2.0 μ m PEEK frit for 102000-P355.



Precolumn Filter

^{**}Other column dimensions and custom-made column dimensions are available.

Ordering Information

How to Order

It's fast and easy to order from the Sepax on-line store at: www.sepax-tech.com

Or contact Sepax Sales Department by

Phone: (302) 366-1101 1-877-SEPAX-US **E** Fax: (302) 366-1151

E-mail:

sales@sepax-tech.com

Mail:

5 Innovation Way, Delaware Technology Park Newark, Delaware 19711 USA



Sepax Technologies, Inc. 5 Innovation Way, Newark, Delaware 19711, USA Tel: (302) 366-1101 Fax: (302) 366-1151 www.sepax-tech.com